




A case series evaluating rates of post-tonsillectomy haemorrhage (PTH) requiring a second surgery in 783 children receiving COBLATION[®] Intracapsular Tonsillectomy (CIT)


Albright JT, Duncan NO, Smerica AM, Edmonds JL. Intra-capsular complete tonsillectomy, a modification of surgical technique to eliminate delayed post-operative bleeding. *Int J Pediatr Otorhinolaryngol.* 2020;128:10970.

Available at: [International Journal of Pediatric Otorhinolaryngology](https://doi.org/10.1016/j.ijot.2020.10970)  

Key points



0 children receiving CIT required a second surgery for PTH over a 12-month period (n=783)




CIT led to significantly less PTH requiring a secondary surgery versus a study of other techniques at the same centre (0 vs 1.41%; p<0.001)

Overview

- A retrospective, multi-surgeon case series evaluating the rate of PTH requiring a second surgery following the introduction of the CIT technique for paediatric tonsillectomy at a single centre in the USA
- A total of 783 children received CIT over a 12-month period
- Both infective and obstructive tonsillectomy indications were included
- The CIT procedure involved attempting a complete tonsillectomy without violating the capsule and injuring the underlying pharyngeal muscles
- PTH rates were compared with a previous study at the same centre (n=7,024) of other tonsillectomy techniques: extracapsular COBLATION, extracapsular COBLATION with partial suture closure, extracapsular diathermy and partial intracapsular tonsillectomy using a microdebrider

Results

- None of the 783 children receiving CIT experienced PTH requiring a second surgery during the 12-month study period
- The rate of PTH requiring a second surgery was significantly lower with CIT than in a previous study of other tonsillectomy techniques at the same centre (Figure)



Conclusions

In a case series of 783 children, no child receiving a COBLATION Intracapsular Tonsillectomy experienced a post-tonsillectomy haemorrhage requiring a second surgery. These findings were significantly improved compared with a previous study at the same centre evaluating other tonsillectomy techniques.