


# Collarless Polished Cemented Stems (CPCS) combined with OXINIUM<sup>◇</sup>/XLPE (VERILAST<sup>◇</sup> Technology) and R3<sup>◇</sup> acetabular cup delivers reliable performance at 10 years


**+ Plus points**



CPCS+OXINIUM/  
XLPE+R3  
Survivorship comparable  
to conventional THR<sup>1</sup>



**10A\***  
ODEP rating for  
CPCS<sup>2</sup>




**10A\***  
ODEP rating for  
R3<sup>2</sup>

**Overview**

- An automated industry report prepared by the Australian Orthopaedic Association National Joint Replacement Registry (AOANJRR) summarises usage and outcomes of CPCS combined with OXINIUM/XLPE bearing and R3 cup compared to all other conventional total hip replacements (THR), based on data collected by the AOANJRR<sup>†</sup>
- CPCS+OXINIUM/XLPE+R3 combination recorded usage between January 2008 and July 2019<sup>1</sup>
  - 2,717 total procedures
  - 2,504 total patients
  - 149 implanting surgeons at 112 centres

**Results**

- Cumulative percent revision of CPCS+OXINIUM/XLPE+R3 combination is similar to conventional THR in the AOANJRR (Figure)<sup>1</sup>



Number at Risk	0 Yr	1 Yr	2 Yrs	3 Yrs	4 Yrs	5 Yrs	6 Yrs	7 Yrs	8 Yrs
CPCS	2717	2243	1837	1482	1174	872	635	425	284
Other Total Conventional Hip	459862	411967	367183	324803	284423	246881	212057	180196	151377

Number at Risk	9 Yrs	10 Yrs	11 Yrs	12 Yrs	13 Yrs	14 Yrs	15 Yrs	16 Yrs	17 Yrs
CPCS	158	65	11	0	0	0	0	0	0
Other Total Conventional Hip	125311	102254	83097	66968	52986	40541	29076	19073	10719

Note: Procedures using metal/metal prostheses with head size larger than 32mm are excluded from the comparator

Figure. Cumulative percent revision of primary total conventional hip replacement by model (all diagnoses)

**Conclusion**

CPCS+OXINIUM/XLPE+R3 demonstrates 10-year survivorship comparable to conventional THR.

<sup>†</sup>AOANJRR is confident in the accuracy of the data included in this report, at the time it was provided. However, it was generated using an automated reporting system and has not been reviewed by the AOANJRR personnel.



## CPCS

Cemented Hip System

**16 years**  
of clinical heritage  
**10A\* ODEP rating**

### Triple-taper design

Trapezoidal stem designed to reduce stresses on the cement mantle,<sup>3</sup> as well as reduce subsidence<sup>4</sup>

### Improved range of motion and stability

Accurate restoration of joint anatomy with standard and high offset options in every size, while the circulo-trapezoidal neck provides an optimised head-to-neck ratio to improve range of motion and stability<sup>5,6</sup>

### Excellent wear performance

Exclusive combination of OXINIUM and highly cross-linked polyethylene (XLPE) has excellent wear performance in laboratory studies, clinical studies and registries<sup>7-10</sup>

### Low levels of taper corrosion

Substantially lower levels of taper corrosion compared to metal femoral heads<sup>11,12</sup>

### Biocompatibility

Contains lower levels of nickel, cobalt and chromium compared to cobalt chromium molybdenum implants<sup>13,14</sup>

### STIKTITE<sup>®</sup> stability

When compared with more traditional porous coatings, STIKTITE coating has greater porosity providing a higher coefficient of friction for an immediate 'scratch-fit' feel and the potential for better initial implant fixation<sup>15,16</sup>

**Improved initial fixation** limits micromotion potentially enhancing bony ingrowth<sup>17</sup>



## OXINIUM<sup>◊</sup>

with XLPE

**15 years**  
of clinical heritage



## R3<sup>◊</sup>

Acetabular System

**15 years**  
of clinical heritage  
**10A\* ODEP rating**

#### References

1. Australian Orthopaedic Association National Joint Replacement Registry (AOANJRR), Automated Industry Report System (AIRS), ID No.1497 for Smith+Nephew, CPCS Total Conventional Hip, (Procedures from 1 September 1999 – 13 August 2019), Accessed 14 August 2019, AOA, Adelaide: 1-19.
2. Orthopaedic Data Evaluation Panel (ODEP). Available at <http://www.odep.org.uk>. Accessed 23 November 2019.
3. Walsh WR, Turner AWL, Bertollo N, Pelletier M, Nielsen G, Gillies RM. Effect of implant size on bone remodeling of cemented hip prosthesis. Poster presented at ORS Congress; March 18-22, 2006; Chicago, IL, USA.
4. McCalden RW, Charron KD, Yuan X, Bourne RB, Naudie DD, MacDonald SJ. Randomised controlled trial comparing early migration of two collarless polished cemented stems using radiostereometric analysis. *J Bone Joint Surg Br*. 2010;92:935-940.
5. Barrack RL. Dislocation after total hip arthroplasty: implant design and orientation. *J Am Acad Orthop Surg*. 2003;11:89-99.
6. Barrack RL, Lavernia C, Ries M, Thornberry R, Tozakoglou E. Virtual reality computer animation of the effect of component position and design on stability after total hip arthroplasty. *Orthop Clin North Am*. 2001;32(4):569-577.
7. Parikh A, Hill P, Pawar V, Sprague J. Long-term simulator wear performance of an advanced bearing technology for THA. Poster 1028 presented at ORS Congress; January 26-29, 2013; San Antonio, TX, USA.
8. Al-Khateeb H, Teeter M, Naudie D, McCalden RW. Oxidised zirconium versus cobalt-chromium in THA: analysis of retrieved acetabular liners. *Orthopaedic Proceedings*. 2016;98-B:(Supp 1)14-14. Abstract only.
9. Jassim S, Patel S, Wardle N, et al. Five-year comparison of wear using oxidised zirconium and cobalt-chrome femoral heads in total hip arthroplasty. *Bone Joint J*. 2015;97:883-889.
10. Australian Orthopaedic Association National Joint Replacement Registry (AOANJRR) Hip, Knee & Shoulder Arthroplasty: 2019 Annual Report. Available at: <https://aoanjrr.sahmri.com/annual-reports-2019>. Accessed 6 January 2020.
11. Li C, Parikh A, Sprague J, Pawar V. Mechanically assisted crevice corrosion on CoCrMo heads after long-term hip simulator wear testing. Poster 0270 presented at ORS Congress March 28-31, 2015; Las Vegas, NV, USA.
12. Cartner J, Aldinger P, Li C, Collins D. Characterization of femoral head taper corrosion features using a 22-year retrieval database. *HSSJ*. 2017;13:35-41.
13. Hunter G, Dickinson J, Herb B, Graham R. Creation of oxidized zirconium orthopaedic implants. *J ATSM Int*. 2005;2:7.
14. ASTM F75-01. Standard specification for cobalt-28 chromium-6 molybdenum alloy castings and casting alloy for surgical implants (UNS R30075), ASTM International, West Conshohocken, PA, 2012.
15. Heiner AD and Brown TD. Frictional coefficients of a new bone ingrowth structure. Poster 1623 presented at ORS Congress; February 11-14, 2007; San Diego, CA, USA.
16. Bourne RB, McCalden RW, Naudie D, Charron KD, Yuan X, Holdsworth DW. The next generation of acetabular shell design and bearing surfaces. *Orthopaedics*. 2008;31(12 Suppl 2):818-826.
17. Pilliar RM. Cementless implant fixation—toward improved reliability. *Orthop Clin N Am*. 2005;36:113-119.