

# Best practices in cementing for better outcomes

Objective: Obtain optimal cement adhesion and penetration for better implant-cement-bone interlock and decreased loosening<sup>1</sup>

## Implant/glove cleanliness:

**Clean/dry implant surfaces increases cement adhesion strength<sup>7</sup>**

1. Use clean/dry gloves to prevent contamination
2. If necessary, quickly clean implant contamination with alcohol to reduce risk (<4 min)

## Bone cleanliness:

**Pulsed lavage increases cement penetration<sup>5</sup> and reduces radiolucent lines<sup>8</sup>**

1. Dry bone after lavage<sup>6</sup>
2. Lavage cleans away blood, debris, bone marrow, and fat, which decreases cement penetration resistance<sup>5</sup>

## Cement penetration:

**3-5mm is the optimal cement mantle<sup>2</sup>**

1. Apply cement equally to tibial component and tibial bone by either spatula or ingerpacking<sup>3,4</sup>
  - 10-16 grams of cement roughly correlates to 3-5mm cement penetration<sup>3</sup>
  - Application via gun injection leads to 5+mm cement mantlen<sup>3</sup>
2. Control bone bleeding (tourniquet, TXA, etc)
  - Bone bleeding can displace cement or create cement voids preventing interlock especially when cement viscosity is low<sup>8</sup>
3. Allow cement to reach dough stage
  - Earlier than dough stage decreases penetration<sup>9</sup>
  - High viscosity cement doughs quickly and the viscosity remains constant during working<sup>10</sup>
4. Apply constant pressure
  - Impact time is short, therefore, cement is not lowing/penetrating bone
  - Extend knee after prosthesis insertion and do not disturb<sup>9</sup>

## Bone survival:

**Cement temperature and quantity affects bone survival**

1. <44°C/111°F for 1 minute = no major changes in bone formation to implant<sup>6</sup>
2. >5mm cement mantle leads to bone necrosis<sup>3,4</sup> due to increased thermal damage



# References

1. Miskovsky C, Whiteside LA, White SE. The cemented unicondylar knee arthroplasty. An in vitro comparison of three cement techniques. *Clin Orthop Relat Res* 1992; 215.
2. Bert JM, McShane M. Is it necessary to cement the tibial stem in cemented total knee arthroplasty? *Clin Orthop Relat Res* 1998; 73.
3. Vanlommel J, Luyckx JP, Labey L, Innocenti B, De Corte R, Bellemans J. Cementing the Tibial Component in Total Knee Arthroplasty: Which Technique is the Best? *JOA* 2011; 26(3).
4. Lutz MJ, Pincus PF, Whitehouse SL, et al. The effect of cement gun and cement syringe use on the tibial cement mantle in total knee arthroplasty. *J Arthroplasty* 2009; 24:46.
5. Clarius M, Seeger JB, Jaeger S, Mohr G, Bitsch RG. The importance of pulsed lavage on interface temperature and ligament tension force in cemented unicompartmental knee arthroplasty. *Clin Biomech* 2012; 27.
6. Eriksson RA, Albrektsson T. The effect of heat on bone regeneration: an experimental study in the rabbit using the bone growth chamber. *J Oral Maxillofac Surg* 1984; 42.
7. Salehi A, Walker B, Choudhury S. Factors influencing the adhesion strength of bone cement to metallic substrates. 56th Annual Meeting of the ORS. Poster No. 2101. 2010.
8. Benjamin JB, Gie GA, Lee JC, Ling RSM, RG Volz. Cementing Technique and the effects of bleeding. *JBJS* 1987; 69.
9. Silverman EJ, Landy DC, Massel DH, Kaimrajh DC, Latta LL, Robinson RP. The effect of viscosity on cement penetration in total knee arthroplasty, an application of the Squeeze Film Effect. *JOA* 2014; 29.
10. Webb JCJ, Spencer RF. The role of polymethylmethacrylate bone cement in modern orthopaedic surgery. *JBJS (Br)* 2007; 89-B.

For detailed product information, including indications for use, contraindications, effects, precautions and warnings, please consult the product's Instructions for Use (IFU) prior to use.

**Supporting healthcare professionals for over 150 years**

---

**Smith & Nephew, Inc.**, 1450 Brooks Road, Memphis, TN 38116, USA  
Phone: 1-901-396-2121, Information: 1-800-821-5700, Orders/inquiries: 1-800-238-7538  
[www.smith-nephew.com](http://www.smith-nephew.com)

®Trademark of Smith & Nephew.

©2016 Smith & Nephew  
05459 V1 02/16