

A SIMPLE SOLUTION

FOR PATIENT POSITIONING



Lateral decubitus
positioning made easy¹

 **smith&nephew**
ACUFEX[®] Access
Positioning System

Supporting healthcare professionals

Customize your lateral decubitus positioning needs without complexity¹

Smith & Nephew's **ACUFEX[®] Access Positioning System** was designed with ease of use in mind.² By providing limb abduction, rotation, flexion and optional lateral traction controls, our dedicated system allows you to customize your positioning solution to meet both technique and repair site needs without complexity.¹ This cost-effective, lightweight and stable solution gives you unimpeded access to the repair site.^{1,2} ACUFEX Access puts the control in the surgeon's hands by providing precise shoulder positioning within the sterile field.

It's Easy.

Simple set up and easy-to-use control points so you can focus on what you do best.¹

It's Modular.

Optional lateral traction makes both standard and advanced positioning configurations possible all while freeing up staff from supporting the operative arm during procedures.

It's Draped.


System drapes allow for surgeon control and intraoperative adjustments within the sterile field eliminating the need for component sterilization.

The **ACUFEX Access Positioning System** includes a cart and is portable with a small storage footprint.



Optional Lateral Traction Boom

Designed to increase access to the repair site. Length and height is adjustable allowing for a range of positions to meet varying patient sizes and surgeon preference.¹



Cable Management

Take the guess work out with a single point of weight attachment and eliminate the need to tie knots with a reusable cable and hook.



Abduction Lever

A simple press of the lever allows surgeons and staff to abduct or adduct the arm within the sterile field.



Internal and External Rotation

Our unique, patent pending rotation control allows surgeons and staff to rotate the arm in one fluid motion within the sterile field. Simply rotate and release to stay in a fixed position.



Universal Clamp

Conveniently attaches to any standard OR table while eliminating the need for adaptors.



Flexion Control

An easy-to-turn hand control is used to flex and extend the arm.¹



Ambidextrous, one size disposable kits. Internal padding provides protection of the neurovascular bundle at the wrist. Unique, patent pending design of the hand plate is anatomically shaped for supporting the hand in a relaxed and natural position.^{2,3}

Ordering information

ACUFEX Access Positioning System

Reference #	Description
72205268	ACUFEX Access Positioning System includes: ACUFEX Access Positioning System ACUFEX Access Cart
72205103	ACUFEX Access Lateral Traction Boom
72205102	ACUFEX Access Positioning Kit
72205104	ACUFEX Access Advanced Positioning Kit



ACUFEX Access Positioning System shown in the advanced setup configuration with patient attached.

To order these instruments, call **+1 800 343 5717** in the U.S. or contact an authorized Smith & Nephew representative.

Sports Medicine
Smith & Nephew, Inc.
150 Minuteman Road
Andover, MA 01810

www.smith-nephew.com
T +978 749 1000
US Customer Service:
+1 800 343 5717

°Trademark of Smith & Nephew.
©2018 Smith & Nephew.
All rights reserved. Printed in USA.
16254 V1 01/19

Supporting healthcare professionals for over 150 years

References

1. Smith & Nephew 2018. Usability Validation Report 15007361 2. Smith & Nephew 2018. Bench Top Validation Report 15007702 3. Loh PY, Yeoh WL, Nakashima H, Muraki S. 2018. *Deformation of the median nerve at different finger postures and wrist angles*. PeerJ 6:e5406