



Combining a small profile with secure cortical fixation

 **smith&nephew**
MINI MAGNUM[®]
Knotless Fixation Implant




MINI MAGNUM[◇] knotless fixation implant


The MINI MAGNUM knotless implant is designed for use in the glenoid and combines a small profile with secure cortical fixation. In addition, the suture ratchet mechanism transports the labral tissue to the glenoid and achieves a strong knotless suture fixation.

Key features

- Independent bone lock, suture tensioning and suture lock provide precision and control
- Incremental and controlled tensioning for precise repair
- Internal automated suture lock eliminates reliance on bone quality for suture lock
- The resulting knotless repair eliminates the need for knots, and their potential for movement-induced knot migration



MINI MAGNUM implant specifications



MINI MAGNUM knotless fixation implant

Reference #	Description
OM-6006	MINI MAGNUM PLUS 2.8MM KNOTLESS IMPLANT SET WITH SPEEDSTITCH [®] MAGNUMWIRE [®] SUTURE CARTRIDGE, WHITE
OM-6007	MINI MAGNUM PLUS 2.8MM KNOTLESS IMPLANT SET WITH SPEEDSTITCH MAGNUMWIRE SUTURE CARTRIDGE, BLUE CO-BRAID
OM-6008	MINI MAGNUM PLUS 2.8MM KNOTLESS IMPLANT SET WITH SPEEDSTITCH MAGNUMWIRE SUTURE CARTRIDGE, BLACK CO-BRAID
OM-9016	DRILL GUIDE AND PATHFINDER [®] SET: COMPATIBLE WITH OM-9080, OM-9310, OM-9316
OM-9018	HIGH VISIBILITY DRILL GUIDE: COMPATIBLE WITH OM-9016, OM-9080, OM-9310, OM-9316
OM-9080	MINI MAGNUM DRILL (3MM)
OM-9332	LABRAFIX [®] METAL INSTRUMENT TRAY
OM-7000	SPEEDSTITCH SUTURING DEVICE HANDLE
OM-8850	SPEEDSTITCH NEEDLE CASSETTE INSERT
OM-8085	SPEEDSTITCH MAGNUMWIRE SUTURE CARTRIDGES, WHITE, 6 PER BOX
OM-8086	SPEEDSTITCH MAGNUMWIRE SUTURE CARTRIDGES, BLUE CO-BRAID, 6 PER BOX
OM-8087	SPEEDSTITCH MAGNUMWIRE SUTURE CARTRIDGES, BLACK CO-BRAID, 6 PER BOX
OM-8088	SPEEDSTITCH MAGNUMWIRE SUTURE CARTRIDGES, MIXED: WHITE, BLUE CO-BRAID AND BLACK CO-BRAID, 6 PER BOX (2 OF EACH COLOR)

ArthroCare Corporation
7000 West William Cannon Drive
Austin, TX 78735
USA

www.smith-nephew.com

Order Entry: 1-800-343-5717
Order Entry Fax: 1-888-994-2782