

 **smith&nephew**

SPEEDLOCK®

Knotless Suture Anchor

For precise, reliable control
of your instability repairs



Providing precise anatomic fixation of soft tissues

Precise, reliable tissue shift and anatomic repair is now achievable with the SPEEDLOCK® Knotless Suture Anchor. This knotless anchor provides independent suture-lock and accurate suture tensioning after anchor insertion.

The SPEEDLOCK Anchor features independent Suture Lock, which allows you to make incremental micro-tensioning suture adjustments for a precise repair. The internal suture locking mechanism eliminates the bone and suture lock variability that may arise from differences in bone quality and arthroscopic knot tying.

- **Inline drill guide** system provides precise, ergonomic, reliable drill placement.
- **Strong fixation:** The SPEEDLOCK Knotless Suture Anchor demonstrates pullout strength of 250N¹.
- **Independent suture tensioning** allows the surgeon to precisely and reliably reduce the tissue to the underlying bone after anchor insertion.
- **Internal suture lock** allows for precise, reproducible suture tensioning, independent of anchor and bone interference.
- **PEEK implant** is radiolucent, MRI-safe, and generally revisable.



1. In 20 lb. bone foam as documented in report 25374-01 Rev. G for the SPEEDLOCK Knotless Suture Anchor.

Ordering Information

SPEEDLOCK Knotless Suture Anchor

| Reference # | Description |
|-------------|----------------------------------|
| OM-7500 | SPEEDLOCK Knotless Suture Anchor |

SPEEDLOCK Inline Drill System

| Reference # | Description |
|-------------|---|
| OM-9816 | 3.5mm Inline Drill Guide |
| OM-9817 | SPEEDLOCK Obturator for Inline Drill Guide |
| OM-9850 | SPEEDLOCK 3.0mm Tapered Twist Drill |
| OM-9860 | SPEEDLOCK 3.0mm Tapered Twist Drill (hard bone) |

SPEEDLOCK Drill Guide

| Reference # | Description |
|-------------|---|
| OM-9715 | SPEEDLOCK Drill Guide, 3.5mm, high visibility |
| OM-9716 | SPEEDLOCK Drill Guide, 3.5mm, low profile |
| OM-9717 | SPEEDLOCK Sharp-tipped Obturator |
| OM-9750 | SPEEDLOCK 3.0mm Tapered Twist Drill |



Supporting healthcare professionals for over 150 years

Manufactured by:
ArthroCare Corporation
7000 West William Cannon Drive
Austin, TX 78735 USA

Distributed by:
Smith & Nephew, Inc.
150 Minuteman Road
Andover, MA 01810 USA

www.smith-nephew.com
T +1 978 749 1000
US Customer Service:
+1 800 343 5717

*Trademark of Smith & Nephew.
©2017 Smith & Nephew. All rights reserved.
Printed in USA. 05616 V2 PN 74027 Rev. B 05/17